

Sling Label/Inspection Instructions

NOTE: Use an indelible ink pen (e.g. Sharpie) to write on the label.

SERIAL #: For factory use only		DATE OF FIRST USE: _____
<p>THIS SLING MUST BE INSPECTED BEFORE USE ON A PATIENT. DISCONTINUE USE IF FABRIC RIPPED OR WORN OUT</p> <p>SCHEDULE INSPECTION EVERY SIX MONTHS TO VERIFY INTEGRITY. WRITE INSPECTION DATES BELOW</p>		
SIX MONTH INSPECTION DATES		
INSPECTION 1	INSPECTION 2	INSPECTION 3
INSPECTION 4	INSPECTION 5	INSPECTION 6

Write date (MM/DD/YY) and name initials when sling is used for the first time.

For each inspection, write date (MM/DD/YY) and inspector's name initials.

Schedule and record sling inspections every 6 months **after** date of first use up to 6 inspections.

It is recommended that the sling be removed from service 6 months after Inspection 6 (3.5 years after Inspection 1).

SLING INSPECTION CHECKLIST:

- ✓ Loops and straps are free of cuts or fraying.
- ✓ Stitching is intact on loops, straps, seams and panels.
- ✓ Fabric has no worn spots, weave separation or visible cuts.

*** REMOVE SLING FROM SERVICE IF IT DOES NOT MEET ONE OR MORE INSPECTION CRITERIA ***