



molift transfer system

molift[®]
designed for life



The best days are those when you feel you have extra energy, either for work or for leisure. The Molift Transfer System is designed to encourage good working positions and relieve your body from strain. The products are easy and comfortable to use, and they function on their own as well as in combination with other products. In the end, their most important function is to get the job done, without putting your back at risk. Invest in your own well being by letting Molift Transfer System help you with the heavy transfers. This will give you extra energy to help more people, help them better, spend quality time with them – or have more energy for yourself!

The philosophy behind the Molift Transfer System

Physical strain may cause irreparable damage. Caregivers are often subject to chronic muscle and skeletal injuries in their daily work with their patients. This results in reduced quality of life and increased absence from work due to these injuries.

The Molift Transfer System consists of simple and user-friendly means of assistance, which reduce load damages while moving patients. Our central concept is that all transfers are done using less body strength, while making it more safe and free of pain for the user. By reducing the friction against the surface and/or providing handles, the transfer can take place easier, more comfortable and effortless for both carer and user. Thereby, the carer can move the user without putting his or her back at risk, while the user can participate, if possible, and hereby retaining his/her physical functions.

The Molift Transfer System is developed in line with modern transfer principles, which are based on:

- the natural movement pattern of the human body
- weight transfer
- the knowledge of the friction points of the body

On the following pages, you will find some examples of common daily transfer situations, which may cause physical damages. It is possible to use other product combinations than those shown.



From sitting to sitting

The user can often take an active part in the movement.

Molift Handyswing Floor and Molift Support

The user may be able to stand on his/her legs.



- Reduces friction
- Provides handles
- Encourages correct working positions
- Allows user to contribute to the transfer



Molift Board and Molift Handytrans

Used on patients who are unable to stand up.



- Reduces friction
- Provides handles
- Encourages correct working positions
- Creates a bridge between wheelchair and bed



From lying to sitting

The user may be able to take an active part in the movement.

Molift Handyswing Seat and Molift Handysafe
Molift Handyswing makes it easy to swing the patient's legs out of bed. Using the Molift Handysafe accessory around the patient's legs, the caregiver can retain a balanced position.



Reduces friction
Encourages the body's natural movement pattern
Creates handles and encourages correct working positions
Allows user to contribute to the transfer
A useful accessory to the Molift Quick Raiser



From lying to lying

Users who are moved from one lying position to another are often incapable of taking any active part in the movement.

Molift Maxitrans and Molift Meditrans
Often used in combination with a draw sheet.
The user is gently transferred from the bed to a stretcher, a shower trolley, etc.



Reduces friction
Padding facilitates a more comfortable and painless transfer
Encourages correct working positions
Easy to put in place and remove

Repositioning

6

Repositioning in bed is one of the tasks, which may cause the greatest risk of damage, and the carer should therefore be extra careful. The combination of the user's natural movement pattern, friction reduction and the carer's active use of weight transfer secures a correct movement situation.



Molift Minitrans and Molift Handysafe

This combination is used when moving the user up in the bed. The user may be able to take some active part in repositioning.



- Easy to put in place and remove
- Reduces friction
- Creates handles
- Encourages a correct working position
- User can contribute to the transfer



Molift Multitrans

Used for moving the user up or sideways in the bed.



- Reduces friction
- Creates handles
- Encourages correct working positions

In and out of a Car

Molift Handyswing Auto

Used to simplify getting into and out of a car.



- Non-slip material underneath prevents it from sliding
- Padding and soft material for extra comfort
- Reasonable alternative to a swing seat
- Neutral colour



Molift Transfer System

| | Product | Primarily used in the following situations: | Example |
|---|--|--|---|
|  | Molift Minitrans Part No. 101060-M | Repositioning in bed | Higher up in the bed |
|  | Molift Meditrans Part No. 101065-M | Transfer from lying to lying | Bed to operating table |
|  | Molift Maxitrans Part No. 101070-M | Transfer from lying to lying | Bed to operating table |
|  | Molift Multitrans Part No. 102040-M | Repositioning in bed | Higher up and sideways in bed |
|  | Molift Handytrans Part No. 102045-M | Repositioning in bed Transfer sitting to sitting | Sideways in bed Bed to wheelchair together with Molift Board |
|  | Molift Handysafe Part No. 104020-M | Supplement to Molift Board, Molift Minitrans and Molift Handyswing Seat. | Bed to chair together with Molift Board |
|  | Molift Support Part No. 104015-M 104014-M | Transfer to a standing position | Up from a chair |
|  | Molift Board Part No. 103060-M 103070-M | Transfer from sitting to sitting | From bed to wheelchair |
|  | Molift Handyswing Seat Part No. 106010-M | Transfer from lying to sitting | To sitting position on the edge of a bed |
|  | Molift Handyswing Floor Part No. 106020-M | Turn in standing position | Bed to chair |
|  | Molift Handyswing Auto Part No. 106015-M | Rotate on a car seat | In and out of a car |
|  | Molift Turn Part No. 101090-M | Change of position in bed | To prevent pressure sores |
|  | Molift One Way Trans Part No. 101050-M | Repositioning in chair | Further into a chair |



Molift - designed for life

At times with increasing demand for efficiency in the healthcare sector and growing focus on patient needs, it is easy to neglect the needs of the carer. The need to give better help and caring avoiding heavy lifting situations, avoiding putting your own back at risk, without being completely exhausted after work and still have some surplus energy when the workday is over.

This is the insight that drives us at Molift. The company which was started 20 years ago by a man who through his own experience saw the need for lifting and moving patients in an effective and comfortable way – for both patient and carer.

We have since then delivered products that compliment the natural movement of the body, that is both simple and intuitive to use. Products designed to give the patient a better life in addition to providing the carer with real support and optimum work enjoyment.

This is the true meaning of; Molift – designed for life

molift[®]

www.molift.com



Molift Ltd.
The Courtyard
Georges Road Corner
Wellington Road North
Stockport, Cheshire SK4 1HT
Tel: 0161 476 2413



Molift Inc.
5008 W. Linebaugh Ave.
Suite 60
Tampa, Florida 33624
Tel: (813) 969-2213

Head Office: Molift AS, Ole Deviks vei 44, 0668 Oslo, Norway, Tel.: +47 4000 1004

For international distributors see www.molift.com