



molift quick raiser

molift[®]
designed for life



The best days are those when you feel you have extra energy-either for work or for leisure. Molift's products are designed to give you this. They are easy and intuitive-both in use and to understand. The design makes it easy to get into narrow and small spaces. In the end, their most important function is to get the job done, without putting your back at risk. Invest in your own well being by letting Molift patient lifters do the heavy work. This will give you extra energy to help more people, help them better, spend quality time with them, or have more energy for yourself.

The philosophy behind Molift Quick Raiser

Molift Quick Raiser is designed for lifting and moving users who can give some assistance during the transfer.

The natural movement pattern is the most energy-saving way to move. Bending forward when standing up is the most natural movement pattern for the body, and the body can use a little help to recognise this. Molift Quick Raiser was designed with this principle in mind. The angled column raises the user in an up-and-forward movement, which imitates, and therefore stimulates, the natural movement pattern. Molift Quick Raiser is the only lift with this unique feature.



Molift Quick Raiser is the best choice for optimum transfers from one sitting position to another, regardless of whether it involves beds, chairs or toilets, standing up or stand-up training. Molift Quick Raiser is available in three different versions.



Molift Easy Raiser secures easy undressing and dressing in toilet situations, with dignity.

Molift Easy Raiser Support

Molift Easy Raiser

The Molift Easy Raiser sling is designed to fit the Molift Quick Raiser and provide quick, comfortable, secure transfers. The sling gives good support for the small of the back and is very easily fitted. The design makes it the perfect sling for toileting. Molift Easy Raiser has an anti-slip surface on the inside to prevent the sling from sliding up. The high back gives an even pressure distribution. The belt around the waist provides extra support and stability. Molift Easy Raiser Support has extra pelvic support for users who need additional help and support to stand. Both slings are available in sizes from XS to XL.

Molift Quick Raiser 1

Ergonomic Design

The Molift Quick Raiser success story began with Molift Quick Raiser 1. It was first introduced in 1996 and was soon rewarded for its understanding of the natural movement patterns of the body and the ergonomic teamwork between nurse and user. One year after its introduction, it was awarded The Norwegian Design Council prize for its pioneering design and functionality.

User-Friendly

Molift Quick Raiser 1 has an extra low base, which makes it easy to get under beds, chairs and most wheelchairs. Even though its legs cannot be spread, they are wide enough to get around a toilet. The small wheels are placed in such a way that the lift is very easy to manoeuvre. Its weight distribution gives stability for safe, secure transfers.

Small and Lightweight

Molift Quick Raiser 1 is small and compact and has very low weight – only 30 kg. It is the smallest active lift on the market, with an incredibly small turning circle, it is equally suitable for narrow places in the home and in hospitals and institutions



The legs of the Molift Quick Raiser 1 can easily get around a toilet

Specification:

Lifting capacity: 160 kg
Leg height: 75 mm
Width: 610 mm
Weight: 30 kg
Part No.: 21011

Molift Quick Raiser 2

Since its first introduction, the Molift Quick Raiser has been the market leader for active lifts, with more than 10,000 units sold in Europe. The angled column, which gives the body a natural movement pattern, has been much appreciated by both nurses and users. It is difficult to match Molift Quick Raiser with regard to efficiency, comfort and dignity, which are valuable characteristics for toileting and all other situations.

Molift Quick Raiser 2 uses the excellent features from Molift Quick Raiser 1, and adds a few more:

Electrical, adjustable legs

If you cannot get under a chair, for example, you can get around it.

Larger Kneepads

Most users feel that larger kneepads are more comfortable. The kneepads can also be turned around if that works better for the user.

Specification:

Lifting capacity:	160 kg
Leg height:	100 mm
Width:	610 mm inner position 1200 mm outer position
Weight:	36 kg
Part No.:	23011



Molift Quick Raiser 2+

Molift Quick Raiser 2+ is the latest addition to the Molift Quick Raiser family. With its increased lifting capacity, it can now lift and transfer a user weighing up to 200 kg.

Specially Designed Lifting Arm

In relation to the two previous models, Molift Quick Raiser 2+ has a specially designed lifting arm, which gives more room for the user.

Low Weight – Large Capacity

Like the others in the family, the Molift Quick Raiser 2+ is surprisingly small and compact and has a low weight of only 37 kg. Like the Molift Quick Raiser 2, the 2+ model has electrical, adjustable legs and has large, self-adjusting kneepads.

All Molift Quick Raisers can be dismantled.



Specification:

Lifting capacity:	200 kg
Leg height:	100 mm
Width:	610 mm inner position 1200 mm outer position
Weight:	37 kg
Part No.:	24011

Molift Quick Raiser for children



Molift Quick Raiser for children

For children and small adults, a height-adjustable footboard can be fitted. It is adjustable to individual needs on both the right and left sides. Maximum user weight is 70 kg.

Accessories:

Sandals for positioning the feet

"Sandals" for both adults and children can be fitted on the footboard of the lift. They stop the feet from sliding or moving during lifts.

Molift Power Pac

The Molift Power Pac system is a battery solution consisting of batteries mounted in an aluminium profile, which is placed in a holder on the column. The battery charger can be wall mounted or a free standing desk top tray is an option, then just plug into a regular wall socket. The charger warranty is two years and up to five years if annual safety checks are done. The battery warranty is one year.

Useful suggestions:

- keep two batteries and have a lift that works all the time;
- always keep a fully charged battery in the charger;

Important information on lifts and movements

Before lifting or moving a user, it is necessary to evaluate the functionality level of the user. You should consider the user's medical and cognitive state, and check for adjustments of his/her surroundings and how motivated the user is for the activity. Based on this, you can decide what technological means of assistance would best support the situation at hand. It is important to take note if the user has functionality fluctuations throughout the day. To make a safe lift and movement, a new evaluation may be needed before each lift. The means of assistance, which worked fine this morning, may not function so well later in the day or in the evening. Molift Quick Raiser is designed for lifting and moving people with the assistance of a qualified professional. Read the user manual carefully for more information.



Molift - designed for life

At times with increasing demand for efficiency in the healthcare sector and growing focus on patient needs, it is easy to neglect the needs of the carer. The need to give better help and caring avoiding heavy lifting situations, avoiding putting your own back at risk, without being completely exhausted after work and still have some surplus energy when the workday is over.

This is the insight that drives us at Molift. The company which was started 20 years ago by a man who through his own experience saw the need for lifting and moving patients in an effective and comfortable way – for both patient and carer.

We have since then delivered products that compliment the natural movement of the body, that is both simple and intuitive to use. Products designed to give the patient a better life in addition to providing the carer with real support and optimum work enjoyment.

This is the true meaning of; Molift – designed for life

molift[®]

www.molift.com

Head Office: Molift Group AS, Ole Deviks vei 44, 0668 Oslo, Norway, Tel.: +47 4000 1004
For international distributors see www.molift.com