



Hovermatt Comparison Table

MODEL	DESCRIPTION	WEIGHT LIMIT	APPLICATION	MAINTENANCE
HM34HS HM39HS HM50HS	HOVERMATT HEAT SEALED 34", 39" AND 50" 	34" 600 LB/ 272 KG	LONG TERM (ICU, FLOOR , BARIATRIC)	WIPE AFTER EACH USE
		39" 850LB/386 KG	MAY REMAIN UNDER PATIENT FOR FULL HOSPITAL STAY	LAUNDRY WASH IF SOILED WITH FLUIDS
		50" 1500 LB/ 680 KG	X-RAY/ MRI ALLOWED	
HM34DC HM39DC	HOVERMATT DOUBLE COATED 34" AND 39" 	34" 600 LB/ 272 KG	HIGH TURNOVER (RADIOLOGY, ER, OR)	WIPE AFTER EACH USE
		39" 850LB/386 KG	MAY ONLY REMAIN UNDER PATIENT FOR MAX. 23 HOURS X-RAY/ MRI ALLOWED	MAY BE WIPED CLEAN IF SOILED WITH FLUIDS